




F90A9 mouse mAb

Catalog No	YP-mAb-12366
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM90A9P FAM90A9
Immunogen	Synthesized peptide derived from human F90A9 AA range: 70-120
Specificity	This antibody detects endogenous levels of F90A9 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F90A9;FAM90A9P FAM90A9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	caution:Could be the product of a pseudogene.,similarity:Belongs to the FAM90 family.,
Background	FAM90A9 belongs to subfamily II of the primate-specific FAM90A gene family, which originated from multiple duplications and rearrangements (Bosch et al., 2007 [PubMed 17684299]). For background information on the FAM90A gene family, as well as information on the evolution of FAM90A genes, see FAM90A1 (MIM 613041).[supplied by OMIM, Oct 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using F90A9 mouse mAb