




MYCT mouse mAb

Catalog No	YP-mAb-12391
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC2A13
Immunogen	Synthesized peptide derived from human MYCT AA range: 352-402
Specificity	This antibody detects endogenous levels of MYCT at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Solute carrier family 2 member 13;Proton myo-inositol cotransporter;Hmit;H(+)-myo-inositol symporter;H(+)-myo-inositol cotransporter;MYCT;SLC2A13
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Predominantly expressed in the brain.
Function	function:H(+)-myo-inositol cotransporter. Can also transport related stereoisomers.,PTM:Glycosylated.,similarity:Belongs to the major facilitator superfamily. Sugar transporter (TC 2.A.1.1) family.,tissue specificity:Predominantly expressed in the brain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MYCT mouse mAb