



CXL10 mouse mAb

Catalog No	YP-mAb-12482
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CXL10
Immunogen	Synthesized peptide derived from mouse CXL10
Specificity	This antibody detects endogenous levels of Mouse CXL10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CXL10;C-X-C motif chemokine 10 (10 kDa interferon gamma-induced protein; Gamma-IP10;IP-10;Small-inducible cytokine B10) [Cleaved into: CXCL10(1-73)]
Calculated Molecular Weight	11kD
Observed Band	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images