



PHGDH mouse mAb

Catalog No	YP-mAb-12503
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	PHGDH PGDH3
Immunogen	Synthesized peptide derived from human PHGDH
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat PHGDH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	3 PGDH;3PGDH;3-PGDH;D-3-phosphoglycerate dehydrogenase;EC:1.1.1.37; PHGDH;PHGDH PGDH3;D-3-phosphoglycerate dehydrogenase (3-PGDH;EC 1.1.1.95)
Observed Band	57 kDa
Calculated Molecular Weight	57 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

