



sesn2 mouse mAb

Catalog No	YP-mAb-12582
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SESN2 SEST2
Immunogen	Synthesized peptide derived from human sesn2
Specificity	This antibody detects endogenous levels of sesn2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.74% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SESN2;Sestrin2;EC:1.11.1.-;Hi95;Hypoxia-induced gene;SESN2 SEST2;Sestrin-2 (Hi95)
Observed Band	54-60 kDa
Calculated Molecular Weight	54 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

