




fxr1 mouse mAb

Catalog No	YP-mAb-12587
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FXR1
Immunogen	Synthesized peptide derived from human fxr1
Specificity	This antibody detects endogenous levels of fxr1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.79% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FMR1 autosomal homolog 1;hFXR1p;RNA-binding protein FXR1;fxr1;Fragile X mental retardation syndrome-related protein 1 (hFXR1p)
Observed Band	68kD
Calculated Molecular Weight	68kD
Cell Pathway	Cytoplasm, cytosol . Adjacent to Z-lines in muscles. .
Tissue Specificity	Expressed in all tissues examined including heart, brain, kidney and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of HeLa and mouse brain lysis, using primary antibody at 1:1000 dilution.

Secondary antibody(catalog#:RS23920) was diluted at 1:10000