



ACSL4 mouse mAb


Catalog No	YP-mAb-12589
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ACSL4 ACS4 FACL4 LACS4
Immunogen	Synthesized peptide derived from human ACSL4
Specificity	This antibody detects endogenous levels of ACSL4 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.81% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACSL4;FACL4;ACS4;ACSL 4;Arachidonate--CoA ligase;ACSL4 ACS4 FACL4 LACS4;Long-chain-fatty-acid--CoA ligase 4 (EC 6.2.1.3) (Long-chain acyl-CoA synthetase 4) (LACS 4)
Observed Band	75kD
Calculated Molecular Weight	75kD
Cell Pathway	Mitochondrion outer membrane ; Single-pass type III membrane protein . Peroxisome membrane ; Single-pass type III membrane protein . Microsome membrane ; Single-pass type III membrane protein . Endoplasmic reticulum membrane ; Single-pass type III membrane protein . Cell membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Regulation of Fatty Acid Metabolism and Inhibition of Colorectal Cancer Progression by Erchen Decoction
Evidence-based Complementary and Alternative Medicine
Linghong Liao, Fei Zhang, Zewei Zhuo, Chengbao Huang, Xiaofang Zhang, Ruifang Liu, Bizhen Gao, Shanshan Ding
WB Mouse colorectal tissue

