




ALS2CR8 Monoclonal Antibody

Catalog No	YP-mAb-12683
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ALS2CR8
Immunogen	The antiserum was produced against synthesized peptide derived from human ALS2CR8. AA range:311-360
Specificity	ALS2CR8 Monoclonal Antibody detects endogenous levels of ALS2CR8 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALS2CR8;Calcium response factor;CARF;NYD SP24;Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 8 protein;Calcium-response factor;Testis development protein NYD-SP24
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus .
Tissue Specificity	Colon,Lymph node,Testis,
Function	function:May be a transcription factor.,
Background	function:May be a transcription factor.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ALS2CR8 Monoclonal Antibody