



APBA2BP Monoclonal Antibody

Catalog No	YP-mAb-12685
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NECAB3
Immunogen	The antiserum was produced against synthesized peptide derived from human NECAB3. AA range:321-370
Specificity	APBA2BP Monoclonal Antibody detects endogenous levels of APBA2BP protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	APBA2BP;EFCBP3;NECAB3;Nek2 interacting protein 1;NIP1;STIP3;SYTIP2;X11L binding protein 51;XB51;N-terminal EF-hand calcium-binding protein 3; Amyloid beta A4 protein-binding family A member 2-binding protein;Nek2-interacting protein 1;Neuronal calcium-binding protein 3;X11L-binding protein 51
Observed Band	44kD
Calculated Molecular Weight	44kD
Cell Pathway	Golgi apparatus .
Tissue Specificity	Strongly expressed in heart and skeletal muscle, moderately in brain and pancreas.
Function	alternative products:Additional isoforms seem to exist. Experimental confirmation may be lacking for some isoforms,function:Inhibits the interaction of APBA2 with beta-amyloid precursor protein (APP), and hence allows formation of beta-amyloid.,PTM:Phosphorylated by NEK2.,similarity:Contains 1 ABM domain.,similarity:Contains 1 EF-hand domain.,subunit:Interacts with the N-terminal domain of APBA2. Interacts with NEK2.,tissue specificity:Strongly expressed in heart and skeletal muscle, moderately in brain and pancreas.,



Background

The protein encoded by this gene interacts with the amino-terminal domain of the neuron-specific X11-like protein (X11L), inhibits the association of X11L with amyloid precursor protein through a non-competitive mechanism, and abolishes the suppression of beta-amyloid production by X11L. This protein, together with X11L, may play an important role in the regulatory system of amyloid precursor protein metabolism and beta-amyloid generation. The protein is phosphorylated by NIMA-related expressed kinase 2, and localizes to the Golgi apparatus. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using APBA2BP Monoclonal Antibody