



Neurexophilin-3 Monoclonal Antibody

Catalog No	YP-mAb-12753
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NXPH3
Immunogen	The antiserum was produced against synthesized peptide derived from human NXPH3. AA range:158-207
Specificity	Neurexophilin-3 Monoclonal Antibody detects endogenous levels of Neurexophilin-3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Neurexophilin-3;NXPH3;KIAA1159;NPH3
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Secreted .
Tissue Specificity	Highest level in brain.
Function	function:May be signaling molecules that resemble neuropeptides. Ligand for alpha-neurexins.,PTM:May be proteolytically processed at the boundary between the N-terminal non-conserved and the central conserved domain in neuron-like cells.,similarity:Belongs to the neurexophilin family.,tissue specificity:Highest level in brain.,
Background	function:May be signaling molecules that resemble neuropeptides. Ligand for alpha-neurexins.,PTM:May be proteolytically processed at the boundary between the N-terminal non-conserved and the central conserved domain in neuron-like cells.,similarity:Belongs to the neurexophilin family.,tissue specificity:Highest level in brain.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using
Neurexophilin-3 Monoclonal Antibody