



# SR-1F Monoclonal Antibody


|                                    |                                                                                                                                                                                                                                                                                                                                                  |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-12813                                                                                                                                                                                                                                                                                                                                     |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                              |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                                                                                  |
| <b>Applications</b>                | WB                                                                                                                                                                                                                                                                                                                                               |
| <b>Gene Name</b>                   | HTR1F                                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen</b>                   | The antiserum was produced against synthesized peptide derived from human 5-HT-1F. AA range:101-150                                                                                                                                                                                                                                              |
| <b>Specificity</b>                 | SR-1F Monoclonal Antibody detects endogenous levels of SR-1F protein.                                                                                                                                                                                                                                                                            |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                                                                          |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                                                                            |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                             |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                                                                                                                                                                                                  |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                          |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                     |
| <b>Synonyms</b>                    | SR-1F;HTR1F;5-hydroxytryptamine receptor 1F;HTR1EL;5-HT-1F;5-HT1F; Serotonin receptor 1F                                                                                                                                                                                                                                                         |
| <b>Observed Band</b>               | 38kD                                                                                                                                                                                                                                                                                                                                             |
| <b>Calculated Molecular Weight</b> | 38kD                                                                                                                                                                                                                                                                                                                                             |
| <b>Cell Pathway</b>                | Cell membrane ; Multi-pass membrane protein .                                                                                                                                                                                                                                                                                                    |
| <b>Tissue Specificity</b>          | Brain,                                                                                                                                                                                                                                                                                                                                           |
| <b>Function</b>                    | function:This is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. The activity of this receptor is mediated by G proteins that inhibit adenylate cyclase activity.,similarity:Belongs to the G-protein coupled receptor 1 family., |
| <b>Background</b>                  | function:This is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. The activity of this receptor is mediated by G proteins that inhibit adenylate cyclase activity.,similarity:Belongs to the G-protein coupled receptor 1 family., |

**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**

Western Blot analysis of various cells using SR-1F Monoclonal Antibody