




GPR107 Monoclonal Antibody

Catalog No	YP-mAb-13283
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR107
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR107. AA range:141-190
Specificity	GPR107 Monoclonal Antibody detects endogenous levels of GPR107 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	G protein coupled receptor 107;GCDRP;KIAA1624;LUSTR1;Protein GPR107; GPR107;Lung seven transmembrane receptor 1
Observed Band	67kD
Calculated Molecular Weight	67kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Golgi apparatus, trans-Golgi network membrane .
Tissue Specificity	Brain,Gastric cancer,Lymph node,Small intestine,
Function	similarity:Belongs to the LU7TM family.,
Background	similarity:Belongs to the LU7TM family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GPR107 Monoclonal Antibody