




GPR108 Monoclonal Antibody

Catalog No	YP-mAb-13284
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR108
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR108. AA range:155-204
Specificity	GPR108 Monoclonal Antibody detects endogenous levels of GPR108 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	G protein coupled receptor 108;GPR108;LUSTR2;Protein GPR108;Lung seven transmembrane receptor 2
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Golgi apparatus, cis-Golgi network membrane ; Multi-pass membrane protein . Golgi apparatus, trans-Golgi network membrane ; Multi-pass membrane protein . Golgi apparatus membrane ; Multi-pass membrane protein . Colocalizes with TLR3, -7, -4, and -9. .
Tissue Specificity	Uterus,
Function	similarity:Belongs to the LU7TM family.,
Background	similarity:Belongs to the LU7TM family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using GPR108 Monoclonal Antibody