




GPR115 Monoclonal Antibody

Catalog No	YP-mAb-13288
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR115
Immunogen	The antiserum was produced against synthesized peptide derived from human GPR115. AA range:641-690
Specificity	GPR115 Monoclonal Antibody detects endogenous levels of GPR115 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPR115;Probable G-protein coupled receptor 115;PGR18;G-protein coupled receptor PGR18
Observed Band	82kD
Calculated Molecular Weight	82kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Testis,
Function	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GPR115 Monoclonal Antibody