



GPR50 Monoclonal Antibody

Catalog No	YP-mAb-13342
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR50
Immunogen	The antiserum was produced against synthesized peptide derived from human MTR1L. AA range:301-350
Specificity	GPR50 Monoclonal Antibody detects endogenous levels of GPR50 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Melatonin related receptor;Mel1c;H9;G protein-coupled receptor 5;G protein coupled receptor 5;GPR50;Melatonin-related receptor;G protein-coupled receptor 50
Observed Band	68kD
Calculated Molecular Weight	68kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Hypothalamus and pituitary.
Function	function:Does not bind melatonin.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Hypothalamus and pituitary.,
Background	This gene product belongs to the G-protein coupled receptor 1 family. Even though this protein shares similarity with the melatonin receptors, it does not bind melatonin, however, it inhibits melatonin receptor 1A function through heterodimerization. Polymorphic variants of this gene have been associated with bipolar affective disorder in women. [provided by RefSeq, Jan 2010],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using GPR50 Monoclonal Antibody