



ALR Monoclonal Antibody

Catalog No	YP-mAb-13741
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GFER
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human GFER. AA range:51-100
Specificity	ALR Monoclonal Antibody detects endogenous levels of ALR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALR_HUMAN;Augmenter of liver regeneration;ERV1;ERV1 homolog;Erv1 like growth factor;FAD-linked sulfhydryl oxidase ALR;GFER;Growth factor augmenter of liver regeneration;Growth factor erv1 like;Hepatic regenerative stimulation substance;Hepatopoietin;Hepatopoietin protein;hERV1;HPO;HPO1;HPO2;HSS;ALR
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	[Isoform 1]: Mitochondrion intermembrane space. Mitochondrion .; [Isoform 2]: Cytoplasm. Secreted.
Tissue Specificity	Ubiquitously expressed. Highest expression in the testis and liver and low expression in the muscle.
Function	catalytic activity:4 R'C(R)SH + O(2) = 2 R'C(R)S-S(R)CR' + 2 H(2)O., cofactor:FAD.,function:FAD-dependent sulfhydryl oxidase that catalyzes disulfide bond formation.,similarity:Contains 1 ERV/ALR sulfhydryl oxidase domain.,subunit:Homodimer.,



Background

The hepatotrophic factor designated augmentor of liver regeneration (ALR) is thought to be one of the factors responsible for the extraordinary regenerative capacity of mammalian liver. It has also been called hepatic regenerative stimulation substance (HSS). The gene resides on chromosome 16 in the interval containing the locus for polycystic kidney disease (PKD1). The putative gene product is 42% similar to the scERV1 protein of yeast. The yeast scERV1 gene had been found to be essential for oxidative phosphorylation, the maintenance of mitochondrial genomes, and the cell division cycle. The human gene is both the structural and functional homolog of the yeast scERV1 gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ALR Monoclonal Antibody