



LSP1 Monoclonal Antibody


Catalog No	YP-mAb-13943
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LSP1
Immunogen	The antiserum was produced against synthesized peptide derived from human LSP1. AA range:104-153
Specificity	LSP1 Monoclonal Antibody detects endogenous levels of LSP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	47 kDa actin-binding protein;Lymphocyte-specific antigen WP34;Lymphocyte-specific protein 1;pp52;WP34;LSP1;52 kDa phosphoprotein
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Cell membrane; Peripheral membrane protein; Cytoplasmic side.
Tissue Specificity	Activated T-lymphocytes.
Function	function:May play a role in mediating neutrophil activation and chemotaxis., PTM:Phosphorylated by casein kinase II, protein kinase C and MAPKAPK2. Phosphorylation by PKC induces translocation from membrane to cytoplasm. Phosphorylation by MAPKAPK2 may regulate neutrophil chemotaxis., subunit:Binds actin.,tissue specificity:Activated T-lymphocytes.,
Background	This gene encodes an intracellular F-actin binding protein. The protein is expressed in lymphocytes, neutrophils, macrophages, and endothelium and may regulate neutrophil motility, adhesion to fibrinogen matrix proteins, and transendothelial migration. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using LSP1 Monoclonal Antibody