




A1BG Monoclonal Antibody

Catalog No	YP-mAb-13971
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	A1BG
Immunogen	Synthesized peptide derived from A1BG . at AA range: 100-180
Specificity	A1BG Monoclonal Antibody detects endogenous levels of A1BG protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A1BG;Alpha 1B-Glycoprotein;A1B;ABG;alpha 1 B glycoprotein;Alpha-1B-glycoprotein;Alpha-1-B glycoprotein
Observed Band	28kD
Calculated Molecular Weight	28kD
Cell Pathway	Secreted.
Tissue Specificity	Plasma.
Function	online information:The Singapore human mutation and polymorphism database, similarity:Contains 5 Ig-like V-type (immunoglobulin-like) domains.,tissue specificity:Plasma.,
Background	The protein encoded by this gene is a plasma glycoprotein of unknown function. The protein shows sequence similarity to the variable regions of some immunoglobulin supergene family member proteins. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using A1BG Monoclonal Antibody