



ADCK1 Monoclonal Antibody


Catalog No	YP-mAb-14643
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ADCK1
Immunogen	The antiserum was produced against synthesized peptide derived from human ADCK1. AA range:251-300
Specificity	ADCK1 Monoclonal Antibody detects endogenous levels of ADCK1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADCK1;Uncharacterized aarF domain-containing protein kinase 1
Observed Band	61kD
Calculated Molecular Weight	61kD
Cell Pathway	Secreted .
Tissue Specificity	Placenta,Skin,
Function	function:The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).,similarity:Belongs to the protein kinase superfamily. ADCK protein kinase family.,similarity:Contains 1 protein kinase domain.,
Background	function:The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).,similarity:Belongs to the protein kinase superfamily. ADCK protein kinase family.,similarity:Contains 1 protein kinase domain.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ADCK1 Monoclonal Antibody