



ADCK4 Monoclonal Antibody

Catalog No	YP-mAb-14645
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ADCK4
Immunogen	The antiserum was produced against synthesized peptide derived from human ADCK4. AA range:31-80
Specificity	ADCK4 Monoclonal Antibody detects endogenous levels of ADCK4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADCK4;Uncharacterized aarF domain-containing protein kinase 4
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	Mitochondrion membrane ; Single-pass membrane protein . Cytoplasm, cytosol . Cell membrane .
Tissue Specificity	Widely expressed, including renal podocytes.
Function	function:The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).,similarity:Belongs to the protein kinase superfamily. ADCK protein kinase family.,similarity:Contains 1 protein kinase domain.,
Background	This gene encodes a protein with two copies of a domain found in protein kinases. The encoded protein has a complete protein kinase catalytic domain, and a truncated domain that contains only the active and binding sites of the protein kinase domain, however, it is not known whether the protein has any kinase activity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ADCK4 Monoclonal Antibody