




NARFL Monoclonal Antibody

Catalog No	YP-mAb-14873
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NARFL
Immunogen	The antiserum was produced against synthesized peptide derived from human NARFL. AA range:387-436
Specificity	NARFL Monoclonal Antibody detects endogenous levels of NARFL protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HPRN;IOP1;LET1L;NARFL;PRN;Protein related to Narf;Cytosolic Fe-S cluster assembly factor NARFL;Iron-only hydrogenase-like protein 1;Nuclear prelamins A recognition factor-like protein
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	CIA complex,
Tissue Specificity	Widely expressed.
Function	function:Regulates HIF1A expression.,similarity:Belongs to the NARF family.,tissue specificity:Widely expressed.,
Background	function:Regulates HIF1A expression.,similarity:Belongs to the NARF family.,tissue specificity:Widely expressed.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using NARFL Monoclonal Antibody