



PDIK1L Monoclonal Antibody

Catalog No	YP-mAb-14899
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PDIK1L
Immunogen	The antiserum was produced against synthesized peptide derived from human PDIK1L. AA range:101-150
Specificity	PDIK1L Monoclonal Antibody detects endogenous levels of PDIK1L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CLIK1L;PDIK1L;RP11 96L14.4;STK35L2;Serine/threonine-protein kinase PDIK1L;PDLIM1-interacting kinase 1-like
Observed Band	39kD
Calculated Molecular Weight	39kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in liver, kidney, pancreas, spleen, thymus and prostate.
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family., similarity:Contains 1 protein kinase domain.,tissue specificity:Expressed in liver, kidney, pancreas, spleen, thymus and prostate.,
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family., similarity:Contains 1 protein kinase domain.,tissue specificity:Expressed in liver, kidney, pancreas, spleen, thymus and prostate.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using PDIK1L Monoclonal Antibody