




PPP4R1L Monoclonal Antibody

Catalog No	YP-mAb-14957
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PPP4R1L
Immunogen	The antiserum was produced against synthesized peptide derived from human PPP4R1L. AA range:191-240
Specificity	PPP4R1L Monoclonal Antibody detects endogenous levels of PPP4R1L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2orf192;PRO185;Protein phosphatase 4 regulatory subunit 1 like;PPP4R1L; Putative serine/threonine-protein phosphatase 4 regulatory subunit 1-like; C20orf192
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	
Tissue Specificity	Fetal liver,Thymus,
Function	function:May be a regulatory subunit of serine/threonine-protein phosphatase 4., similarity:Contains 4 HEAT repeats.,
Background	function:May be a regulatory subunit of serine/threonine-protein phosphatase 4., similarity:Contains 4 HEAT repeats.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using PPP4R1L Monoclonal Antibody