



# Hox-C6 Monoclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                     |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-15768                                                                                                                                                                                                                                                                                        |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                 |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                                     |
| <b>Applications</b>                | WB                                                                                                                                                                                                                                                                                                  |
| <b>Gene Name</b>                   | HOXC6                                                                                                                                                                                                                                                                                               |
| <b>Immunogen</b>                   | The antiserum was produced against synthesized peptide derived from human HOXC6. AA range:186-235                                                                                                                                                                                                   |
| <b>Specificity</b>                 | Hox-C6 Monoclonal Antibody detects endogenous levels of Hox-C6 protein.                                                                                                                                                                                                                             |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                                                             |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                                                                                                                                               |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                                                                                                                                                     |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                             |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                        |
| <b>Synonyms</b>                    | CP25;HHO.C8;Homeo box 3C;Homeo box C6;Homeo box C8 protein; Homeobox C6;Homeobox protein CP25;Homeobox protein HHO.C8;Homeobox protein Hox-3C;Homeobox protein Hox-C6;HOX3;HOX3C;Hoxc6; HXC6_HUMAN;Hox-C6                                                                                           |
| <b>Observed Band</b>               | 26kD                                                                                                                                                                                                                                                                                                |
| <b>Calculated Molecular Weight</b> | 26kD                                                                                                                                                                                                                                                                                                |
| <b>Cell Pathway</b>                | Nucleus.                                                                                                                                                                                                                                                                                            |
| <b>Tissue Specificity</b>          | Lung,                                                                                                                                                                                                                                                                                               |
| <b>Function</b>                    | function:Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.,similarity:Belongs to the Antp homeobox family., similarity:Contains 1 homeobox DNA-binding domain.,     |
| <b>Background</b>                  | This gene belongs to the homeobox family, members of which encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located |



on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene, HOXC6, is one of several HOXC genes located in a cluster on chromosome 12. Three genes, HOXC5, HOXC4 and HOXC6, share a 5' non-coding exon. Transcripts may include the shared exon spliced to the gene-specific exons, or they may include only the gene-specific exons. Alternatively spliced transcript variants encoding different isoforms have been identified for HOXC6. Transcript variant two includes the shared exon, and transcript variant one includes only gene-specific exons. [provided by RefSeq, Jul 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images



Western Blot analysis of various cells using Hox-C6 Monoclonal Antibody