



GIPR Monoclonal Antibody


Catalog No	YP-mAb-15919
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GIPR
Immunogen	The antiserum was produced against synthesized peptide derived from human GIPR. AA range:93-142
Specificity	GIPR Monoclonal Antibody detects endogenous levels of GIPR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Glucose-dependent insulintropic polypeptide receptor;GIP-R;GIP R;GIPR; Gastric inhibitory polypeptide receptor
Observed Band	53kD
Calculated Molecular Weight	53kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Colon,G-protein coupled receptors,Pancreas,
Function	function:This is a receptor for GIP. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.,similarity:Belongs to the G-protein coupled receptor 2 family.,
Background	This gene encodes a G-protein coupled receptor for gastric inhibitory polypeptide (GIP), which was originally identified as an activity in gut extracts that inhibited gastric acid secretion and gastrin release, but subsequently was demonstrated to stimulate insulin release in the presence of elevated glucose. Mice lacking this gene exhibit higher blood glucose levels with impaired initial insulin response after oral glucose load. Defect in this gene thus may contribute to the pathogenesis of diabetes. [provided by RefSeq, Oct 2011],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using GIPR Monoclonal Antibody