



IGFBP7 Mouse mAb

Catalog No	YP-mAb-15926
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-282 of human IGFBP7
Specificity	
Formulation	
Source	Monoclonal, Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IBP7;IBP-7;IGF-binding protein 7;IGFBP-7;Insulin-like growth factor-binding protein 7;IGFBP7;MAC25;PSF;IGFBP-rP1;MAC25 protein;PGI2-stimulating factor;Prostacyclin-stimulating factor;Tumor-derived adhesion factor;TAF;AGM;FSTL2;RAMSVPS;IGFBP-7v;IGFBPRP1
Observed Band	35kDa
Calculated Molecular Weight	29kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	This gene encodes a member of the insulin-like growth factor (IGF)-binding protein (IGFBP) family. IGFBPs bind IGFs with high affinity, and regulate IGF availability in body fluids and tissues and modulate IGF binding to its receptors. This protein binds IGF-I and IGF-II with relatively low affinity, and belongs to a subfamily of low-affinity IGFBPs. It also stimulates prostacyclin production and cell adhesion. Alternatively spliced transcript variants encoding different isoforms have been described for this gene, and one variant has been



associated with retinal arterial macroaneurysm (PMID:21835307).


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of various lysates using IGFBP7 Mouse pAb (A13023) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 90s.