



CGREF1 Monoclonal Antibody

Catalog No	YP-mAb-16026
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CGREF1
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CGREF1. AA range:1-50
Specificity	CGREF1 Monoclonal Antibody detects endogenous levels of CGREF1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGRE1;CGR11;Cell growth regulator with EF hand domain protein 1;CGREF1;Cell growth regulatory gene 11 protein;Hydrophobestin
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	Secreted .
Tissue Specificity	Brain,Fetal brain,Testis,Thalamus,
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Able to inhibit growth in several cell lines.,induction:By p53.,similarity:Contains 2 EF-hand domains.,
Background	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Able to inhibit growth in several cell lines.,induction:By p53.,similarity:Contains 2 EF-hand domains.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CGREF1 Monoclonal Antibody