




HRP-3 Monoclonal Antibody

Catalog No	YP-mAb-16040
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HDGFRP3
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human HDGFRP3. AA range:131-180
Specificity	HRP-3 Monoclonal Antibody detects endogenous levels of HRP-3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HDGF2;HRP-3;HDGF-2;CGI-142;HDGFRP3;Hepatoma-derived growth factor-related protein 3;Hepatoma-derived growth factor 2
Observed Band	23kD
Calculated Molecular Weight	23kD
Cell Pathway	Nucleus .
Tissue Specificity	Detected in testis, heart, spinal cord and brain.
Function	function:Enhances DNA synthesis and may play a role in cell proliferation., similarity:Belongs to the HDGF family.,similarity:Contains 1 PWWP domain., tissue specificity:Detected in testis, heart, spinal cord and brain.,
Background	function:Enhances DNA synthesis and may play a role in cell proliferation., similarity:Belongs to the HDGF family.,similarity:Contains 1 PWWP domain., tissue specificity:Detected in testis, heart, spinal cord and brain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using HRP-3 Monoclonal Antibody