




IL-22 Monoclonal Antibody

Catalog No	YP-mAb-16075
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IL22
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human IL22. AA range:121-170
Specificity	IL-22 Monoclonal Antibody detects endogenous levels of IL-22 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IL-22;If2b1;Iltif;IL-22a;ILTIFa;IL22;Interleukin-22;ZCYTO18;Cytokine Zcyto18;IL-10-related T-cell-derived-inducible factor;IL-TIF
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	Secreted.
Tissue Specificity	Blood,PCR rescued clones,
Function	function:Cytokine that contributes to the inflammatory response in vivo., similarity:Belongs to the IL-10 family.,
Background	function:Cytokine that contributes to the inflammatory response in vivo., similarity:Belongs to the IL-10 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using IL-22
Monoclonal Antibody