



CTACK Monoclonal Antibody

Catalog No	YP-mAb-16077
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCL27
Immunogen	Synthesized peptide derived from C-C motif chemokine 27 at AA range: 51-100
Specificity	CTACK Monoclonal Antibody detects endogenous levels of CTACK protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALP;C C motif chemokine 27;CC chemokine ILC;CCL27;CTACK;CTAK;ESKINE;IL 11 R alpha locus chemokine;ILC;PESKY;SCYA27;Skinkine;Small inducible cytokine A27;C-C motif chemokine 27;Cutaneous T-cell-attracting chemokine;IL-11 R-alpha-locus chemokine;Small-inducible cytokine A27
Observed Band	12kD
Calculated Molecular Weight	12kD
Cell Pathway	Secreted .
Tissue Specificity	Testis, thymus, placenta, ovary and skin.
Function	function:Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10.,online information:CCL27 entry,similarity:Belongs to the intercrine beta (chemokine CC) family.,tissue specificity:Testis, thymus, placenta, ovary and skin.,
Background	This gene is one of several CC cytokine genes clustered on the p-arm of chromosome 9. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The protein encoded by this gene is



chemotactic for skin-associated memory T lymphocytes. This cytokine may also play a role in mediating homing of lymphocytes to cutaneous sites. It specifically binds to chemokine receptor 10 (CCR10). Studies of a similar murine protein indicate that these protein-receptor interactions have a pivotal role in T cell-mediated skin inflammation. [provided by RefSeq, Sep 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CTACK Monoclonal Antibody