




GIMAP2 Monoclonal Antibody

Catalog No	YP-mAb-16163
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GIMAP2
Immunogen	The antiserum was produced against synthesized peptide derived from human GIMAP2. AA range:201-250
Specificity	GIMAP2 Monoclonal Antibody detects endogenous levels of GIMAP2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GIMAP2;GTPase IMAP family member 2;GTPase;IMAP family member 2; HIMAP2;IMAP2;Immunity associated protein 2;Immunity-associated protein 2
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Lipid droplet .
Tissue Specificity	Detected in T-cells.
Function	similarity:Belongs to the IAN GTP-binding protein family.,
Background	This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using GIMAP2 Monoclonal Antibody