




GRTP1 Monoclonal Antibody

Catalog No	YP-mAb-16166
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GRTP1
Immunogen	The antiserum was produced against synthesized peptide derived from human GRTP1. AA range:71-120
Specificity	GRTP1 Monoclonal Antibody detects endogenous levels of GRTP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ22474;Growth hormone-regulated TBC protein 1;Grtp1;GRTP1_HUMAN;MGC138328;MGC13833;OTTHUMP18752;TBC1 domain family member 6;TBC1 domain family;member 6;TBC1D6
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	intracellular,endomembrane system,
Tissue Specificity	Brain,Cervix carcinoma,Kidney epithelium,Ovary,Pancreas,Uterus,
Function	function:May act as a GTPase-activating protein for Rab family protein(s)., similarity:Contains 1 Rab-GAP TBC domain.,
Background	function:May act as a GTPase-activating protein for Rab family protein(s)., similarity:Contains 1 Rab-GAP TBC domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using GRTP1 Monoclonal Antibody