




RAB24 Mouse mAb

Catalog No	YP-mAb-16191
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ras related protein Rab 24;Ras-related protein Rab-24;Rab 24;RAB24
Observed Band	23kDa
Calculated Molecular Weight	23kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	RAB24 is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking (Olkkonen et al., 1993 [PubMed 8126105]).
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunoprecipitation analysis of 200 μ g extracts of MCF-7 cells, using 3 μ g RAB24 antibody (A8394). Western blot was performed from the immunoprecipitate using RAB24 antibody (A8394) at a dilution of 1:1000.