



Ran BP-6 Monoclonal Antibody


Catalog No	YP-mAb-16223
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RANBP6
Immunogen	The antiserum was produced against synthesized peptide derived from human RANBP6. AA range:800-849
Specificity	Ran BP-6 Monoclonal Antibody detects endogenous levels of Ran BP-6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C631B19;FLJ11112;FLJ54311;OTTHUMP2139;OTTHUMP214;RAN binding protein 6;Ran GTP binding protein;RANBP 6;Ran BP-6;RANBP6;Ran-binding protein 6
Observed Band	125kD
Calculated Molecular Weight	125kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Endometrium,Kidney,Placenta,
Function	function:May function in nuclear protein import as nuclear transport receptor., similarity:Belongs to the importin beta family.,similarity:Contains 7 HEAT repeats.,
Background	function:May function in nuclear protein import as nuclear transport receptor., similarity:Belongs to the importin beta family.,similarity:Contains 7 HEAT repeats.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Ran BP-6 Monoclonal Antibody