




RHOBTB3 Monoclonal Antibody

Catalog No	YP-mAb-16246
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RHOBTB3
Immunogen	The antiserum was produced against synthesized peptide derived from human RHOBTB3. AA range:335-384
Specificity	RHOBTB3 Monoclonal Antibody detects endogenous levels of RHOBTB3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0878;RHOBTB3;Rho-related BTB domain-containing protein 3
Observed Band	70kD
Calculated Molecular Weight	70kD
Cell Pathway	Golgi apparatus .
Tissue Specificity	Ubiquitous. Highly expressed in neural and cardiac tissues, pancreas, placenta and testis.
Function	similarity:Contains 2 BTB (POZ) domains.,tissue specificity:Ubiquitous. Highly expressed in neural and cardiac tissues, pancreas, placenta and testis.,
Background	RHOBTB3 is a member of the evolutionarily conserved RHOBTB subfamily of Rho GTPases. For background information on RHOBTBs, see RHOBTB1 (MIM 607351).[supplied by OMIM, Apr 2004],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RHOBTB3 Monoclonal Antibody