




RRP22 Monoclonal Antibody

Catalog No	YP-mAb-16248
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RASL10A
Immunogen	The antiserum was produced against synthesized peptide derived from human RASL10A. AA range:107-156
Specificity	RRP22 Monoclonal Antibody detects endogenous levels of RRP22 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ras like protein family member 1A;Ras like protein RRP22;RAS like;family 1; member A;RAS related on chromosome 22;Ras related protein on chromosome 22;Ras related protein on chromosome 22 homolog;Ras-like protein family member 1A;Ras-like protein RRP22;Ras-related protein on chromosome 22; RASL 1A;Rasl1a;RASL1A protein;Rasl1a RAS like;RRP 22;RRP22; RSLAA_HUMAN;RASL10A;Ras-like protein family member 10A
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side . Nucleus, nucleolus . May cycle in and out of the nucleolus in a GTP-dependent manner.
Tissue Specificity	Expression appears to be strictly limited to the central nervous system.
Function	similarity:Belongs to the small GTPase superfamily. Ras family.,
Background	similarity:Belongs to the small GTPase superfamily. Ras family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RRP22 Monoclonal Antibody