




Cdc50B Monoclonal Antibody

Catalog No	YP-mAb-16687
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM30B
Immunogen	The antiserum was produced against synthesized peptide derived from human TMEM30B. AA range:91-140
Specificity	Cdc50B Monoclonal Antibody detects endogenous levels of Cdc50B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	91311B11Rik;CC5B_HUMAN;CDC5;S. cerevisiae;homolog of;B;CDC5B;Cell cycle control protein 5B;MGC126775;TMEM3B;Transmembrane protein 3B;Cdc50B;TMEM30B;Cell cycle control protein 50B;Transmembrane protein 30B
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Heart,Lung,
Function	similarity:Belongs to the CDC50/LEM3 family.,
Background	similarity:Belongs to the CDC50/LEM3 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using Cdc50B Monoclonal Antibody