



MCM7 Monoclonal Antibody

Catalog No	YP-mAb-16795
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MCM7
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human MCM7. AA range:491-540
Specificity	MCM7 Monoclonal Antibody detects endogenous levels of MCM7 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	P1.1 MCM3;MCM2;EC:3.6.4.12;CDC47 homolog;CDC47;MCM7;DNA replication licensing factor MCM7;P1.1-MCM3
Observed Band	81kD
Calculated Molecular Weight	81kD
Cell Pathway	Nucleus . Chromosome . Associated with chromatin before the formation of nuclei and detaches from it as DNA replication progresses. .
Tissue Specificity	Brain,Cervix,Chronic myeloid leukemia cell,Lung,
Function	function:Acts as a factor that allows the DNA to undergo a single round of replication per cell cycle. Required for DNA replication and cell proliferation. Required for S-phase checkpoint activation upon UV-induced damage., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the MCM family.,similarity:Contains 1 MCM domain., subunit:Interacts with the ATR-ATRIP complex and with RAD17. Interacts with TIPIN.,



Background

The protein encoded by this gene is one of the highly conserved mini-chromosome maintenance proteins (MCM) that are essential for the initiation of eukaryotic genome replication. The hexameric protein complex formed by the MCM proteins is a key component of the pre-replication complex (pre_RC) and may be involved in the formation of replication forks and in the recruitment of other DNA replication related proteins. The MCM complex consisting of this protein and MCM2, 4 and 6 proteins possesses DNA helicase activity, and may act as a DNA unwinding enzyme. Cyclin D1-dependent kinase, CDK4, is found to associate with this protein, and may regulate the binding of this protein with the tumorsuppressor protein RB1/RB. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MCM7 Monoclonal Antibody