



Connexin-26 Monoclonal Antibody

Catalog No	YP-mAb-17007
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GJB2
Immunogen	The antiserum was produced against synthesized peptide derived from human Connexin-26. AA range:45-94
Specificity	Connexin-26 Monoclonal Antibody detects endogenous levels of Connexin-26 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GJB2;Connexin 26;CX26;DFNA3;DFNA3A;Connexin-26;Gap junction beta-2 protein
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Cell junction, gap junction . Colocalizes with GJB4 at gap junction plaques in the cochlea. .
Tissue Specificity	Blood,Colon,
Function	disease:Defects in GJB2 are a cause of keratitis-ichthyosis-deafness syndrome (KID syndrome) [MIM:148210]; an autosomal dominant form of ectodermal dysplasia. Ectodermal dysplasias (EDs) constitute a heterogeneous group of developmental disorders affecting tissues of ectodermal origin. EDs are characterized by abnormal development of two or more ectodermal structures such as hair, teeth, nails and sweat glands, with or without any additional clinical sign. Each combination of clinical features represents a different type of ectodermal dysplasia. KID syndrome is characterized by the association of hyperkeratotic skin lesions with vascularizing keratitis and profound sensorineural hearing loss. Clinical features include deafness, ichthyosis,



photobia, absent or decreased eyebrows, sparse or absent scalp hair, decreased sweating and dysplastic finger and toenails.,disease:Defects in GJB2 ar

Background

This gene encodes a member of the gap junction protein family. The gap junctions were first characterized by electron microscopy as regionally specialized structures on plasma membranes of contacting adherent cells. These structures were shown to consist of cell-to-cell channels that facilitate the transfer of ions and small molecules between cells. The gap junction proteins, also known as connexins, purified from fractions of enriched gap junctions from different tissues differ. According to sequence similarities at the nucleotide and amino acid levels, the gap junction proteins are divided into two categories, alpha and beta. Mutations in this gene are responsible for as much as 50% of pre-lingual, recessive deafness. [provided by RefSeq, Oct 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Connexin-26 Monoclonal Antibody