



M-cadherin Monoclonal Antibody

Catalog No	YP-mAb-17048
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CDH15
Immunogen	The antiserum was produced against synthesized peptide derived from human CDH15. AA range:81-130
Specificity	M-cadherin Monoclonal Antibody detects endogenous levels of M-cadherin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAD15_HUMAN;Cadherin 14;Cadherin 15;type 1;M cadherin (myotubule); Cadherin 3;Cadherin-14;Cadherin-15;Cadherin14;Cadherin15;Cadherin3;CCAD; CDH 14;CDH 15;CDH 3;CDH14;Cdh15;CDH3;CDHM;M-cadherin;MCAD;MRD3; Muscle cadherin;Myotubule cadherin
Observed Band	89kD
Calculated Molecular Weight	89kD
Cell Pathway	Cell membrane; Single-pass type I membrane protein.
Tissue Specificity	Expressed in the brain and cerebellum.
Function	disease:A chromosomal aberration involving CDH15 and KIRREL3 is found in a patient with severe mental retardation and dysmorphic facial features. Translocation t(11;16)(q24.2;q24).,disease:Defects in CDH15 are the cause of mental retardation autosomal dominant type 3 (MRD3) [MIM:612580]. Mental retardation is characterized by significantly sub-average general intellectual functioning associated with impairments in adaptive behavior and manifested during the developmental period.,function:Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic



manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. M-cadherin is part of the myogenic program and may provide a trigger for terminal muscle differentiation.,similarity:Contains 5 cadherin domains.,tissue specificity:Expressed in the brai

Background

This gene is a member of the cadherin superfamily of genes, encoding calcium-dependent intercellular adhesion glycoproteins. Cadherins consist of an extracellular domain containing 5 cadherin domains, a transmembrane region, and a conserved cytoplasmic domain. Transcripts from this particular cadherin are expressed in myoblasts and upregulated in myotubule-forming cells. The protein is thought to be essential for the control of morphogenetic processes, specifically myogenesis, and may provide a trigger for terminal muscle cell differentiation. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using M-cadherin Monoclonal Antibody