



HMGCR Mouse mAb

Catalog No	YP-mAb-17137
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB,IP
Gene Name	HMGCR
Immunogen	A synthetic peptide of human HMGCR
Specificity	Cardiovascular
Host	Mouse
Purification	Affinity Purified
Dilution	WB: 1/500-1/1000 IP: 1/20
Concentration	Unconjugated
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	Purification:Affinity Purified
Calculated Molecular Weight	Isotype:IgG
Cell Pathway	
Tissue Specificity	
Function	Buffer System:50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Background	Background:Transmembrane glycoprotein that is the rate-limiting enzyme in cholesterol biosynthesis as well as in the biosynthesis of nonsterol isoprenoids that are essential for normal cell function including ubiquinone and geranylgeranyl proteins.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

