



PPAR alpha mouse mAb

Catalog No	YP-mAb-17767
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PPARA
Immunogen	The antiserum was produced against synthesized peptide derived from human PPAR-alpha. AA range:6-55
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1G1E10;hPPAR;NR1C1;Peroxisome proliferator-activated receptor alpha;PPAR;PPARA; PPARA;PPAR-alpha;Nuclear receptor subfamily 1 group C member 1
Observed Band	Calculated MW: 52 kDa; Observed MW: 52 kDa
Calculated Molecular Weight	Calculated MW: 52 kDa; Observed MW: 52 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Ligand-activated transcription factor.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

