



NLRP3 mouse mAb

Catalog No	YP-mAb-17806
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NLRP3
Immunogen	A synthesized peptide derived from human NALP3
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NLRP3;FCU;MWS;FCAS;Cias1;Mmig1;Pypaf1;AII/AVP;AGTAVPRL
Observed Band	Calculated MW: 118 kDa; Observed MW: 118 kDa
Calculated Molecular Weight	Calculated MW: 118 kDa; Observed MW: 118 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	As the sensor component of the NLRP3 inflammasome, plays a crucial role in innate immunity and inflammation.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

