



MTAP Mouse mAb

Catalog No	YP-mAb-18431
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MTAP MSAP
Immunogen	Synthesized peptide derived from human MTAP
Specificity	This antibody detects endogenous levels of MTAP at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BDMF;MSAP;DMSFH;LGMBF;DMSMFH;c86fus;HEL-249;MTAP;MTAP MSAP; S-methyl-5'-thioadenosine phosphorylase (5'-methylthioadenosine phosphorylase) (MTA phosphorylase) (MTAP) (MTAPase)
Observed Band	29kDa
Calculated Molecular Weight	31kDa
Cell Pathway	Cytoplasm. Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	Catalyzes the reversible phosphorylation of S-methyl-5'-thioadenosine (MTA) to adenine and 5-methylthioribose-1-phosphate. Involved in the breakdown of MTA, a major by-product of polyamine biosynthesis. Responsible for the first step in the methionine salvage pathway after MTA has been generated from S-adenosylmethionine. Has broad substrate specificity with 6-aminopurine nucleosides as preferred substrates.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images