



LRRC47 Mouse mAb

Catalog No	YP-mAb-18435
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LRRC47 KIAA1185
Immunogen	Synthesized peptide derived from human LRRC47
Specificity	This antibody detects endogenous levels of LRRC47 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Leucine-rich repeat-containing protein 47;KIAA1185;LRRC47 KIAA1185
Observed Band	63 kDa
Calculated Molecular Weight	583 aa, 63 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
