



Notch1 Mouse mAb

Catalog No	YP-mAb-18439
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 2456-2555 of human Notch1 (P46531).
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500 - 1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hN1;AOS5;TAN1;AOVD1;Notch1
Observed Band	120kDa/300kDa
Calculated Molecular Weight	273kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	This gene encodes a member of the NOTCH family of proteins. Members of this Type I transmembrane protein family share structural characteristics including an extracellular domain consisting of multiple epidermal growth factor-like (EGF) repeats, and an intracellular domain consisting of multiple different domain types. Notch signaling is an evolutionarily conserved intercellular signaling pathway that regulates interactions between physically adjacent cells through binding of Notch family receptors to their cognate ligands. The encoded preproprotein is proteolytically processed in the trans-Golgi network to generate two polypeptide chains that heterodimerize to form the mature cell-surface receptor. This receptor plays a role in the development of numerous cell and



tissue types. Mutations in this gene are associated with aortic valve disease, Adams-Oliver syndrome, T cell acute lymphoblastic leukemia, chronic lymphocytic leukemia, and head and neck squamous cell carcinoma.


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunoprecipitation analysis of 200 μ g extracts of 293T cells using 3 μ g Notch1 antibody (A19090). Western blot was performed from the immunoprecipitate using Notch1 antibody (A19090) at a dilution of 1:1000