



Band 3/AE1 Mouse mAb

Catalog No	YP-mAb-18527
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC4A1 AE1 DI EPB3
Immunogen	Synthesized peptide derived from human Band 3/AE1
Specificity	This antibody detects endogenous levels of Band 3/AE1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AE 1;AE1;Anion exchange protein 1;Anion exchanger 1;B3AT_HUMAN;Band 3 anion transport protein;BND3;CD233;DI;Diego blood group;EMPB3;EPB3;Erythrocyte membrane protein band 3;Erythroid anion exchange protein;FR;Froese blood group;RTA1A;SLC4A1;Solute carrier family 4 anion exchanger member 1;Solute carrier family 4 member 1;SW;Swann blood group;Waldner blood group;WD;WD1;WR;Wright blood group;SLC4A1 AE1 DI EPB3;Band 3 anion transport protein (Anion exchange protein 1) (AE 1) (Anion exchanger 1) (Solute carrier family 4 member 1) (CD antigen CD233)
Observed Band	100kD
Calculated Molecular Weight	12 kDa; 12kD(Calculated).
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Basolateral cell membrane ; Multi-pass membrane protein . Detected in the erythrocyte cell membrane and on the basolateral membrane of alpha-intercalated cells in the collecting duct in the kidney. .
Tissue Specificity	Detected in erythrocytes (at protein level). ; [Isoform 2]: Expressed in kidney (at protein level).

**Function**

Functions both as a transporter that mediates electroneutral anion exchange across the cell membrane and as a structural protein. Major integral membrane glycoprotein of the erythrocyte membrane; required for normal flexibility and stability of the erythrocyte membrane and for normal erythrocyte shape via the interactions of its cytoplasmic domain with cytoskeletal proteins, glycolytic enzymes, and hemoglobin. Functions as a transporter that mediates the 1:1 exchange of inorganic anions across the erythrocyte membrane. Mediates chloride-bicarbonate exchange in the kidney, and is required for normal acidification of the urine.

Background**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images