



CAPZB Mouse mAb

Catalog No	YP-mAb-18532
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CAPZB
Immunogen	Synthesized peptide derived from human CAPZB
Specificity	This antibody detects endogenous levels of CAPZB at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAPB;CAPZ;CAPPB;CAPZB;F-actin-capping protein subunit beta (CapZ beta)
Observed Band	31kDa
Calculated Molecular Weight	31kDa
Cell Pathway	Cytoplasm, cytoskeleton . Cytoplasm, myofibril, sarcomere .
Tissue Specificity	
Function	F-actin-capping proteins bind in a Ca(2+)-independent manner to the fast growing ends of actin filaments (barbed end) thereby blocking the exchange of subunits at these ends. Unlike other capping proteins (such as gelsolin and severin), these proteins do not sever actin filaments. Plays a role in the regulation of cell morphology and cytoskeletal organization.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images