



Loricrin Mouse mAb

Catalog No	YP-mAb-18571
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LOR LRN
Immunogen	Synthesized peptide derived from human Loricrin
Specificity	This antibody detects endogenous levels of Loricrin at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LOR;LORICRIN;LOR LRN
Observed Band	26kDa
Calculated Molecular Weight	26kDa
Cell Pathway	Cytoplasm. Nucleus, nucleoplasm.
Tissue Specificity	
Function	Major keratinocyte cell envelope protein.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

