




FSTL3 Mouse mAb

Catalog No	YP-mAb-18590
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-263 of human FSTL3 (NP_005851.1)
Specificity	
Formulation	
Source	
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLRG;FSTL3;Follistatin-related protein 3;Follistatin-related gene protein; Follistatin-like protein 3;FSRP
Observed Band	28-33 kDa
Calculated Molecular Weight	263 aa, 28 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Follistatin-like 3 is a secreted glycoprotein of the follistatin-module-protein family. It may have a role in leukemogenesis.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of extracts of various cell lines, using FSTL3 antibody (A10279) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 60s.